GC EUROPE N.V.
Interleuvenlaan 13
B-3001 Leuven

MANUFACTURED by
GC EUROPE N.V.
Interleuvenlaan 13, B-3001 Leuven, Belgium

DISTRIBUTED by
GC EUROPE N.V.
Research Park, Interleuvenlaan 13, B-3001 Leuven,
Belgium

GC AMERICA INC.
3737 West 127th Street, Alsip, ILL. 60803 U.S.A.

PRINTED IN E.U.

GC Initial INmetalbond

Prior to use, carefully read the instructions

GC Initial INmetalbond

INmetalbond is used as a thin layer between the alloy and the first Opaque layer.

The INmetalbond blocs the escaping metal oxides and neutralizes differences in the expansion coefficient. The bonder does not increase the bonding strength as such but allows a wider span of CTE compatibility.

Recommended indications

- Can be used on all PFM alloys, precious as well as non-precious.

Directions for Use

- Metal framework needs to be prepared according manufacturer's instructions.
- Oxidation firing again according manufacturer's instructions.
- Stir paste before using.
- Apply the Bonder in a thin layer but mask the metal completely.
IMPORTANT:
To change the consistency of the INmetalbond paste, use the "GC Initial Paste Opaque Thinner". Only use it in very small quantities.
To avoid drying out, close the cap after using.
- The fired Bonder should have a yellowish, slightly shiny surface.
NOTE:
The colour of the fired Bonder can vary depending on the composition of the alloys.
When using NPA (non-precious alloys) with the INmetalbond, the final temperature of the Opaque Washbake does not need to be increased by 20°C anymore.
e.g.: 1st Opaque Firing on NPA:
GC Initial MC = 940°C
GC Initial LF = 830°C

FIRING INSTRUCTION

| Preheating Temp. | Drying Time | T° increase | Vacuum | Final Temp. | Holding Time |
|------------------|-------------|-------------|--------|-------------|-------------------|
| 550°C | 6min | 80°C/min | Yes | 980°C | 1min no vacuum |

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